


ARRANGEMENT OF RING ELECTRODE CONTACTS ON FLEXIBLE OR PUSHABLE TUBES FOR RECORDING AND STIMULATION. PLEASE CONTACT US WITH DESIGN REQUESTS!

Example of a Percutaneous Electrode with 4 ring-shaped contacts.



General Information:

- Learn more about our technological competences on page 29 and page 34.
- For more information about connection options and materials see pages 30ff.

PRODUCT CONFIGURATION

Number of contacts, length and pitch are arbitrary.
Tube diameter starting from 0.5 mm.

CONNECTION OPTIONS

Optimal interconnection for this product will be specified according to the individual layout and application.

Please see page 31 for further information!